

TILDAS Dual Laser CO₂ Isotope Analyzer for $\delta^{13}\text{C}$, $\delta^{18}\text{O}$ and $\Delta^{17}\text{O}$

Direct Spectroscopic Measurement of $\Delta^{17}\text{O}$ and $\delta^{13}\text{C}$ in CO₂ for Geochemistry and Atmospheric Chemistry

Features


- Direct measurement of CO₂ isotopes in air with no chemical processing or separation
- Repeatability < 0.013 ‰ in $\Delta^{17}\text{O}$ and < 0.01 ‰ in $\delta^{13}\text{C}$ for a 30-minute measurement including a balanced working reference
- Small sample size (~0.25 μmol CO₂ or 15 mL of air) for discrete sample model
- Suitable for samples derived from carbonate via acid digestion
- Automated valve control capable of custom sample scheduling, backgrounds, and calibrations
- 10 Hz time response and 1-sec δ precisions < 0.2 ‰ enable eddy covariance studies for continuous flow model

TILDAS TECHNOLOGY

Aerodyne instruments use **tunable infrared laser direct absorption spectroscopy (TILDAS)** at mid-IR wavelengths to probe molecules at their strongest “fingerprint” transition frequencies. We further enhance sensitivity by employing a patented multi-pass broad-band absorption cell that provides optical path lengths up to 400 meters. Direct absorption spectroscopy allows for fast (<1 sec) absolute trace gas concentrations without need for elaborate calibration procedures. Moreover, TILDAS instruments are relatively free of measurement interference from other molecular species, enabling extremely specific detection.

Related Instruments

- Single laser isotope monitor for $\delta^{13}\text{C}$ and $\delta^{18}\text{O}$ of CO₂
- Single laser isotope monitor for $\Delta^{17}\text{O}$ of CO₂
- Dual laser isotope monitor for CO₂ ($\delta^{13}\text{C}$, $\delta^{18}\text{O}$) and H₂O ($\delta^{18}\text{O}$, δD)



Rugged, field-ready instruments

Direct absorption spectroscopy allows for highly specific and accurate gas detection

Mid-IR detection enables maximum measurement sensitivity


Applications

- Analysis of CO₂ samples derived from carbonate
- Determination of atmospheric sources, sinks, and transport through CO₂ isotopic ratios
- Carbon capture and sequestration monitoring
- Breath analysis

Advantages over IRMS

- Direct measurement of ¹⁷O-CO₂, which is not possible by IRMS
- High precision (< 0.013 ‰ in $\Delta^{17}\text{O}$ and < 0.01 ‰ in $\delta^{13}\text{C}$)
- Lower cost
- Faster measurements (e.g. 30 min)

Powerful software provides easy, flexible instrument control, and real-time result.



TILDAS Dual Laser CO₂ Isotope Analyzer for $\delta^{13}\text{C}$, $\delta^{18}\text{O}$ and $\Delta^{17}\text{O}$

Specifications

Discrete Samples – High Precision

	CO ₂	$\delta^{13}\text{CO}_2$	$\delta\text{CO}^{18}\text{O}$	$\delta\text{CO}^{17}\text{O}$	$\Delta^{17}\text{O}$
1 Sample (0.25 $\mu\text{mol CO}_2$, 3 min)	0.02 ppm	0.03 ‰	0.03 ‰	0.04 ‰	0.04 ‰
10 Samples (2.5 $\mu\text{mol CO}_2$, 30 min)	0.01 ppm	0.01 ‰	0.01 ‰	0.013 ‰	0.013 ‰

Note: These measurements alternate the sample gas with a working reference with a similar mixing ratio, and the time to do so is included in the quoted measurement time.

Continuous Air Measurement – High Precision

	CO ₂	$\delta^{13}\text{CO}_2$	$\delta\text{CO}^{18}\text{O}$	$\delta\text{CO}^{17}\text{O}$
2 min	0.02 ppm	0.03 ‰	0.03 ‰	0.04 ‰
20 min	0.01 ppm	0.01 ‰	0.01 ‰	0.013 ‰

Note: These measurements are normalized to a working reference with a mixing ratio similar to the sample. The flow rate is 0.6 slpm.

Continuous Air Measurement – High Speed

	CO ₂	$\delta^{13}\text{CO}_2$	$\delta\text{CO}^{18}\text{O}$	$\delta\text{CO}^{17}\text{O}$
0.1 second	0.15 ppm	0.3 ‰	0.3 ‰	0.4 ‰
1 second	0.05 ppm	0.1 ‰	0.1 ‰	0.13 ‰
60 second	0.015 ppm	0.03 ‰	0.03 ‰	0.04 ‰

Note: These measurements are not referenced to a working reference. This configuration supports 10 Hz eddy covariance measurements with a modest pump (120 lpm) and a flow rate of 6 slpm.

Dynamic Range

0 – 1000 ppm CO₂


Operating Conditions

-20 to +40 °C, 0 to 20 slpm

REFERENCES

- Vincent J. Hare*, Christoph Dyroff, David D. Nelson, and Drake A. Yarian, High-Precision Triple Oxygen Isotope Analysis of Carbon Dioxide by Tunable Infrared Laser Absorption Spectroscopy, *American Chemical Society, Anal. Chem.* 2022, 94, 46, 16023–16032, October 24, 2022.
- Nathan Perdue, Zachary Sharp, David Nelson, Rick Wehr, Christoph Dyroff, A rapid high-precision analytical method for triple oxygen isotope analysis of CO₂ gas using tunable infrared laser direct absorption spectroscopy, *Rapid Communication in Mass Spectrometry*, 03 September 2022.
- Nitzan Yanay, Rapid and precise measurement of carbonate clumped isotopes using laser spectroscopy, *Science Advances*, 26 Oct 2022, Vol 8, Issue 43
- Steur, P. M., Scheeren, H. A., Nelson, D. D., McManus, J. B., and Meijer, H. A. J., Simultaneous measurement of $\delta^{13}\text{C}$, $\delta^{18}\text{O}$ and $\delta^{17}\text{O}$ of atmospheric CO₂ – Performance assessment of a dual-laser absorption spectrometer. *Atmospheric Measurement Techniques*, 14(6), pp. 4279–4304, 2021.
- McManus, J. Barry, David D. Nelson, and Mark S. Zahniser. Design and performance of a dual-laser instrument for multiple isotopologues of carbon dioxide and water, *Optics Express*, 23(5), pp. 6569–6586, 2015.

Experimental Spectrum Acquired at 1 Hz



Installation

Benchtop system

Flushing the optics with CO₂-free gas is recommended

Instrument Components

Core instrument, thermoelectric chiller, vacuum pump (customer specified), inlet sampling system (customizable), keyboard, mouse, monitor

Data Outputs

RS-232, USB drive, ethernet

Size, Weight, Power

Dimensions: 440 mm x 660 mm x 6U (267mm)
(W x D x H)

Weight: 35 kg (core instrument) + 15 kg (chiller) + pump weight

Electrical Power: 250 W, 120/240 V, 50/60 Hz
(without pump)