

TILDAS Dual Laser CO₂ Isotope Analyzer for $\delta^{13}\text{C}$, $\delta^{18}\text{O}$ and $\Delta^{17}\text{O}$

Direct Spectroscopic Measurement of $\Delta^{17}\text{O}$ and $\delta^{13}\text{C}$ in CO₂ for Geochemistry and Atmospheric Chemistry

Features

- Direct measurement of CO₂ isotopes in air with no chemical processing or separation
- Repeatability < 0.013 ‰ in $\Delta^{17}\text{O}$ and < 0.01 ‰ in $\delta^{13}\text{C}$ for a 30-minute measurement including a balanced working reference
- Small sample size (~0.25 μmol CO₂ or 15 mL of air) for discrete sample model
- Suitable for samples derived from carbonate via acid digestion
- Automated valve control capable of custom sample scheduling, backgrounds, and calibrations
- 10 Hz time response and 1-sec δ precisions < 0.2 ‰ enable eddy covariance studies for continuous flow model

TILDAS TECHNOLOGY

Aerodyne instruments use **tunable infrared laser direct absorption spectroscopy (TILDAS)** at mid-IR wavelengths to probe molecules at their strongest “fingerprint” transition frequencies. We further enhance sensitivity by employing a patented multi-pass broad-band absorption cell that provides optical path lengths up to 400 meters. Direct absorption spectroscopy allows for fast (<1 sec) absolute trace gas concentrations without need for elaborate calibration procedures. Moreover, TILDAS instruments are relatively free of measurement interference from other molecular species, enabling extremely specific detection.

Related Instruments

- Single laser isotope monitor for $\delta^{13}\text{C}$ and $\delta^{18}\text{O}$ of CO₂
- Single laser isotope monitor for $\Delta^{17}\text{O}$ of CO₂
- Dual laser isotope monitor for CO₂ ($\delta^{13}\text{C}$, $\delta^{18}\text{O}$) and H₂O ($\delta^{18}\text{O}$, δD)



Rugged, field-ready instruments

Direct absorption spectroscopy allows for highly specific and accurate gas detection

Mid-IR detection enables maximum measurement sensitivity


Applications

- Analysis of CO₂ samples derived from carbonate
- Determination of atmospheric sources, sinks, and transport through CO₂ isotopic ratios
- Carbon capture and sequestration monitoring
- Breath analysis

Advantages over IRMS

- Direct measurement of ^{17}O -CO₂, which is not possible by IRMS
- High precision (< 0.013 ‰ in $\Delta^{17}\text{O}$ and < 0.01 ‰ in $\delta^{13}\text{C}$)
- Lower cost
- Faster measurements (e.g. 30 min)

Powerful software provides easy, flexible instrument control, and real-time result.



TILDAS Dual Laser CO₂ Isotope Analyzer for $\delta^{13}\text{C}$, $\delta^{18}\text{O}$ and $\Delta^{17}\text{O}$

Specifications

Discrete Samples – High Precision

	CO ₂	$\delta^{13}\text{CO}_2$	$\delta\text{CO}^{18}\text{O}$	$\delta\text{CO}^{17}\text{O}$	$\Delta^{17}\text{O}$
1 Sample (0.25 $\mu\text{mol CO}_2$, 3 min)	0.02 ppm	0.03 ‰	0.03 ‰	0.04 ‰	0.04 ‰
10 Samples (2.5 $\mu\text{mol CO}_2$, 30 min)	0.01 ppm	0.01 ‰	0.01 ‰	0.013 ‰	0.013 ‰

Note: These measurements alternate the sample gas with a working reference with a similar mixing ratio, and the time to do so is included in the quoted measurement time.

Continuous Air Measurement – High Precision

	CO ₂	$\delta^{13}\text{CO}_2$	$\delta\text{CO}^{18}\text{O}$	$\delta\text{CO}^{17}\text{O}$
2 min	0.02 ppm	0.03 ‰	0.03 ‰	0.04 ‰
20 min	0.01 ppm	0.01 ‰	0.01 ‰	0.013 ‰

Note: These measurements are normalized to a working reference with a mixing ratio similar to the sample. The flow rate is 0.6 slpm.

Continuous Air Measurement – High Speed

	CO ₂	$\delta^{13}\text{CO}_2$	$\delta\text{CO}^{18}\text{O}$	$\delta\text{CO}^{17}\text{O}$
0.1 second	0.15 ppm	0.3 ‰	0.3 ‰	0.4 ‰
1 second	0.05 ppm	0.1 ‰	0.1 ‰	0.13 ‰
60 second	0.015 ppm	0.03 ‰	0.03 ‰	0.04 ‰

Note: These measurements are not referenced to a working reference. This configuration supports 10 Hz eddy covariance measurements with a modest pump (120 lpm) and a flow rate of 6 slpm.

Dynamic Range

0 – 1000 ppm CO₂

Operating Conditions

-20 to +40 °C, 0 to 20 slpm

REFERENCES

Vincent J. Hare*, Christoph Dyroff, David D. Nelson, and Drake A. Yarian, High-Precision Triple Oxygen Isotope Analysis of Carbon Dioxide by Tunable Infrared Laser Absorption Spectroscopy, American Chemical Society, Anal. Chem. 2022, 94, 46, 16023–16032, October 24, 2022.


Nathan Perdue, Zachary Sharp, David Nelson, Rick Wehr, Christoph Dyroff, A rapid high-precision analytical method for triple oxygen isotope analysis of CO₂ gas using tunable infrared laser direct absorption spectroscopy, Rapid Communication in Mass Spectrometry, 03 September 2022.

Nitzan Yanay, Rapid and precise measurement of carbonate clumped isotopes using laser spectroscopy, Science Advances, 26 Oct 2022, Vol 8, Issue 43

Steur, P. M., Scheeren, H. A., Nelson, D. D., McManus, J. B., and Meijer, H. A. J., Simultaneous measurement of $\delta^{13}\text{C}$, $\delta^{18}\text{O}$ and $\delta^{17}\text{O}$ of atmospheric CO₂ – Performance assessment of a dual-laser absorption spectrometer. Atmospheric Measurement Techniques, 14(6), pp. 4279–4304, 2021.

McManus, J. Barry, David D. Nelson, and Mark S. Zahniser. Design and performance of a dual-laser instrument for multiple isotopologues of carbon dioxide and water, Optics Express, 23(5), pp. 6569–6586, 2015.

Experimental Spectrum Acquired at 1 Hz



Installation

Benchtop system

Flushing the optics with CO₂-free gas is recommended

Instrument Components

Core instrument, thermoelectric chiller, vacuum pump (customer specified), inlet sampling system (customizable), keyboard, mouse, monitor

Data Outputs

RS-232, USB drive, ethernet

Size, Weight, Power

Dimensions: 440 mm x 660 mm x 6U (267mm)
(W x D x H)

Weight: 35 kg (core instrument) + 15 kg (chiller) + pump weight

Electrical Power: 250 W, 120/240 V, 50/60 Hz
(without pump)